

Tryptophan metabolism Analysis Results

Client	HEM pharma
Company Institution	HEM pharma
Report No.	HEM-ARR-204
Reported Date	2026-03-31
Number of Samples	34

경기 수원시 영통구 창룡대로 256번길 50
광교플렉스데시앙 716호 (16229)

Tel. 031-214-9421
Fax. 070-8680-9311

CONTENTS

	Tryptophan		00 Sample information
			01 Sample preparation
			02 Instrument information
			03 Metabolite analysis results
	Appendix		01 Sample preparation and Data analysis

00. Sample information

Sample type	Feces
Number of Samples	34
Analysis target	Tryptophan metabolism
Analytical instrument	LC-MS/MS

» LC-MS/MS : Liquid Chromatography Triple Quadrupole Mass Spectrometry

01. Sample preparation

No.	Sample ID	Weight (mg)	50% Methanol (mL)
1	FD-599	48.75	0.975
2	FD-602	49.30	0.986
3	FD-605	48.64	0.973
4	FD-606	49.08	0.982
5	FD-607	49.49	0.990
6	FD-610	49.85	0.997
7	FD-616	48.36	0.967
8	FD-622	49.52	0.990
9	FD-624	48.35	0.967
10	FD-629	48.63	0.973
11	FD-642	49.58	0.992
12	FD-646	48.92	0.978
13	FD-648	48.61	0.972
14	FD-649	49.49	0.990
15	FD-706	49.40	0.988
16	FD-737	49.61	0.992
17	FD-762	48.49	0.970
18	FD-765	48.93	0.979
19	FD-775	48.55	0.971
20	FD-790	49.86	0.997
21	FD-791	49.08	0.982
22	FD-803	49.70	0.994
23	FD-868	48.91	0.978
24	FD-916	48.80	0.976
25	FD-920	50.00	1.000
26	FD-921	48.76	0.975
27	FD-926	48.95	0.979
28	FD-945	48.23	0.965
29	FD-952	49.49	0.990
30	FD-1016	48.23	0.965
31	FD-1033	49.58	0.992
32	FD-1046	49.15	0.983
33	FD-1078	49.12	0.982
34	FD-1079	48.64	0.973

» Weight (mg) : 분석에 사용한 시료의 양

» Methanol (mL) : 분석에 사용한 50% 메탄올의 양

02. Instrument information

Analytical instrument	Shimadzu LCMS-8050
Column information	ACQUITY UPLC HSS T3 Column (1.8 μ m, 2.1 mm X 150 mm)
Run time	12 min / sample
Column temperature	40 °C
Mobile phase	(A) Water with 0.1% (v/v) Formic acid (B) Acetonitrile with 0.1% (v/v) Formic acid
Flow rate	0.3 mL / min
Injection volume	1 μ L




Figure 1. Gradient system used in the LC-MS/MS method

» (Fig. 1) 분석 시간(Run time)에 따른 이동상(Mobile phase) 변화

03. Metabolite analysis results


Sample information		Concentration (μmol/g)																	
No.	Sample ID	APAP	AP	AA	DA	GABA	GLU	HA	IND	IAM	IAA	ILA	IPA	KYNA	NA	PA	TA	TRP	XA
1	FD-599	0.000	0.003	0.000	0.000	0.000	3.165	0.000	0.532	0.000	0.029	0.000	0.002	0.000	0.083	0.002	0.125	0.091	0.001
2	FD-602	0.000	0.000	0.000	0.001	0.653	0.434	0.000	0.278	0.000	0.008	0.008	0.010	0.000	0.061	0.000	0.001	0.527	0.005
3	FD-605	0.000	0.000	0.000	0.000	0.000	5.363	0.000	0.569	0.000	0.008	0.002	0.001	0.000	0.072	0.012	0.004	0.946	0.004
4	FD-606	0.000	0.003	0.000	0.000	0.000	3.044	0.000	0.524	0.000	0.002	0.000	0.000	0.001	0.018	0.001	0.233	0.000	0.001
5	FD-607	0.245	0.000	0.000	0.000	0.000	8.588	0.000	2.076	0.000	0.058	0.000	0.001	0.004	0.215	0.000	0.000	1.295	0.008
6	FD-610	0.000	0.001	0.000	0.000	0.181	3.545	0.000	0.520	0.000	0.003	0.000	0.002	0.000	0.100	0.000	0.000	0.619	0.004
7	FD-616	0.000	0.001	0.000	0.000	0.074	7.819	0.000	0.777	0.000	0.060	0.000	0.004	0.002	0.221	0.000	0.000	0.652	0.004
8	FD-622	0.000	0.001	0.000	0.000	0.000	10.419	0.000	0.830	0.000	0.119	0.004	0.007	0.003	0.331	0.003	0.000	0.275	0.005
9	FD-624	0.000	0.001	0.000	0.000	0.147	7.544	0.000	1.914	0.000	0.023	0.000	0.007	0.003	0.155	0.000	0.000	0.987	0.007
10	FD-629	0.000	0.001	0.000	0.000	0.298	1.237	0.000	0.750	0.000	0.000	0.000	0.004	0.001	0.090	0.000	0.002	0.467	0.003
11	FD-642	0.000	0.000	0.000	0.000	0.896	6.061	0.000	1.246	0.000	0.012	0.005	0.002	0.000	0.161	0.000	0.214	1.497	0.009
12	FD-646	0.000	0.001	0.000	0.000	3.312	0.466	0.000	0.226	0.000	0.003	0.000	0.004	0.000	0.071	0.000	0.163	1.716	0.008
13	FD-648	0.002	0.000	0.000	0.000	0.000	1.358	0.000	0.327	0.000	0.019	0.000	0.001	0.002	0.029	0.000	0.000	0.165	0.002
14	FD-649	0.000	0.002	0.000	0.000	0.000	3.567	0.000	0.893	0.000	0.000	0.000	0.001	0.000	0.065	0.000	0.000	0.312	0.002
15	FD-706	0.000	0.002	0.000	0.000	0.000	6.202	0.000	1.272	0.000	0.016	0.005	0.001	0.000	0.112	0.000	0.006	0.298	0.011
16	FD-737	0.000	0.001	0.000	0.000	0.000	5.463	0.000	0.559	0.000	0.007	0.001	0.008	0.003	0.285	0.000	0.000	0.445	0.005
17	FD-762	0.000	0.000	0.000	0.000	0.311	3.911	0.000	0.469	0.000	0.063	0.000	0.002	0.001	0.091	0.000	0.000	0.634	0.003
18	FD-765	0.000	0.001	0.000	0.000	0.000	3.669	0.000	1.628	0.000	0.012	0.000	0.007	0.003	0.020	0.000	0.002	1.318	0.007
19	FD-775	0.000	0.000	0.000	0.000	0.324	6.015	0.000	0.815	0.000	0.054	0.000	0.003	0.000	0.133	0.000	0.005	0.863	0.006
20	FD-790	0.000	0.000	0.000	0.000	0.044	2.895	0.000	0.244	0.000	0.000	0.000	0.001	0.000	0.074	0.000	0.053	0.497	0.003
21	FD-791	0.000	0.000	0.000	0.000	0.000	7.103	0.000	1.371	0.000	0.022	0.003	0.004	0.001	0.231	0.000	0.000	0.716	0.004
22	FD-803	0.000	0.000	0.000	0.000	0.060	10.412	0.000	0.576	0.000	0.015	0.003	0.002	0.002	0.091	0.003	0.000	0.723	0.005
23	FD-868	0.000	0.002	0.000	0.001	0.681	4.404	0.000	1.078	0.000	0.002	0.001	0.004	0.001	0.135	0.000	0.008	0.297	0.002
24	FD-916	0.502	0.000	0.000	0.000	0.483	5.261	0.000	0.596	0.000	0.013	0.000	0.004	0.003	0.194	0.000	0.003	1.471	0.011
25	FD-920	0.000	0.000	0.004	0.000	0.000	2.993	0.000	0.352	0.000	0.009	0.000	0.003	0.000	0.009	0.000	0.000	0.804	0.004
26	FD-921	0.000	0.001	0.000	0.000	0.059	4.338	0.000	0.798	0.000	0.000	0.000	0.003	0.001	0.084	0.000	0.000	0.872	0.005
27	FD-926	0.000	0.000	0.000	0.000	0.000	2.652	0.000	0.501	0.000	0.104	0.000	0.002	0.002	0.047	0.000	0.000	0.375	0.002
28	FD-945	0.000	0.000	0.000	0.001	0.093	4.412	0.000	0.259	0.000	0.002	0.000	0.003	0.001	0.116	0.000	0.000	0.514	0.004

29	FD-952	0.000	0.000	0.000	0.000	0.000	4.350	0.000	0.578	0.000	0.000	0.007	0.001	0.000	0.037	0.002	0.000	0.462	0.003
30	FD-1016	0.000	0.000	0.000	0.000	0.000	3.225	0.000	0.282	0.000	0.000	0.000	0.000	0.002	0.026	0.000	0.000	0.427	0.003
31	FD-1033	0.000	0.000	0.000	0.000	0.000	5.653	0.000	1.059	0.000	0.005	0.001	0.002	0.005	0.088	0.000	0.000	0.490	0.004
32	FD-1046	0.000	0.001	0.000	0.000	0.077	8.208	0.000	0.818	0.000	0.015	0.002	0.002	0.003	0.174	0.000	0.000	0.539	0.004
33	FD-1078	0.003	0.000	0.000	0.000	1.337	7.532	0.000	0.713	0.000	0.003	0.002	0.009	0.003	0.116	0.000	0.000	1.006	0.008
34	FD-1079	0.000	0.000	0.000	0.000	0.000	5.473	0.000	0.182	0.000	0.000	0.000	0.008	0.001	0.121	0.000	0.000	0.705	0.005

» **Abbreviations** : APAP, Acetaminophen; AP, Aminophenol; AA, Anthranilic acid; DA, Dopamine; GABA, gamma-Aminobutyric acid; GLU, Glutamic acid; HA, Histamine; IND, Indole; IAM, Indole-acetamide; IAA, Indole-acetic acid; ILA, Indole-lactic acid; IPA, Indolepropionic acid; KYNA, Kynurenic acid; NA, Nicotinic acid; PA, Picolinic acid; TA, Tryptamine; TRP, Tryptophan; XA, Xanthurenic acid

» **Concentration (μmol/g)** : 분변 g 당 검출된 각 담즙산의 농도(μmol)

01. Sample preparation and Data analysis



- **Sample** : 시료는 입고 즉시 시료의 변질을 최소화 하기 위해 실험 전까지 deep-freezer(-80 °C)에 보관합니다.
- **Metabolite extraction** : 혈청 볼륨에 4배에 해당하는 100% 메탄올을 첨가하여 충분히 볼텍싱 후 다음 원심분리를 합니다. 상층액 필터 후, 96 well plate에 옮깁니다
- **Instrument analysis** : Instrument information 에 따라 LC-MS/MS분석을 수행합니다.
- **Data processing** : 혈청에서 검출된 대사체에 대한 calibration curve를 구축하고, 상관계수 $r > 0.999$ 의 값을 확인하며 해당 curve를 기반으로 분변 내 각 대사체에 대한 정량 값을 산출합니다.

Disclaimer